




# BVT TECHNOLOGIES, a.s.

ELECTROCHEMICAL SENSORS AND DEVICES

*New Conductivity Sensor CC4*




**Autumn 2022**



**Dear customers,**

**We are pleased to inform you about  
our new Conductivity Sensor the  
CC4.**



*„Electrochemical sensors and  
biosensors can be effectively used  
for biochemical activity measure-  
ment.“*

*Dr. Jan Krejci, CEO*

*Front of Sensor*



*Back of Sensor*




**CC4**

## *Uses of the CC4 Conductivity Sensor*

The CC4 conductivity sensor has many uses, and we will list some here to give a good understanding of what it can be used for.

- **Checking the quality of distillation.**
- **Control of water quality in labs (distilled water must have low conductivity)**
- **Control of dialysing solution during dialysis in Kidney Failure Treatment.**
- **Checking of waste waters from treatment plants, i.e. checking of Salt Content**
- **Checking of water from old mines.**
- **Checking of water conductivity in Hydroelectric Dams; if the conductivity will be high, there is a risk of damage to the Dams mechanical parts and the structure.**

## *CC4 Dimensions*



**A— 25.1mm**

**B—26.3mm**

**C—19.3mm**

**D—7.26mm**


## *CC4 Basic Structure*

The Basic Structure of the CC4 is as shown Below. The Structure consists of two Platinum Working Electrodes (Pt WE) and an integrated Thermistor on the opposite side covered by a Dielectric Layer. This Sensor is then integrated with a specialised Cell to create the CC4.


conductivity measurement




temperature measurement



## *Example CC4 Measurement Setup*



## Measurement Software Setup



Recording in PSTrace 5.7, impedance spectroscopy measurement method


The measured electroimpedance spectroscopy (EIS) records from PSTrace 5.7 (dependence of log resistance [ $\Omega$ ] vs. log frequency [Hz]) were transferred to MSExcel.

In the measured records, the area of a constant signal at a frequency value of 20 kHz was determined (corresponds to the measured point of the log frequency of 4.2924).

The subtracted value of the resistance [ $\Omega$ ] corresponding to this frequency was subsequently converted to an electrical conductance of  $1/R$  [S].

The result is the calibration curve of the sensors, as the dependence of the measured electrical conductance of the sensor  $1/R$  [S] in four calibration solutions KCl 1; 0,1; 0,01 and 0,001 mol/l on the specific conductivity of KCl [S/m] standards [1, 2]. From the slope of the calibration curves of the sensors, the cell constant of the sensor is determined, allowing the measurement of conductivity [S/m] in unknown solutions.

## Electrical conductance of the sensor (at a frequency of 20 kHz) on the specific conductivity of KCl calibration solutions



Sample of measured EIS recording of sensor 1 in PStace

Evaluation in Excel




## Results

Electrical conductance of the sensor (at a frequency of 20 kHz) on the specific conductivity of KCl calibration solutions

Electric conductance 1/R [S]	Concentration KCl [mol/l]	Conductivity of KCl standards [S/m]
2.3623E-05	0,001	0,0147
1.5710E-04	0,01	0,1413
1.2468E-03	0,1	1,2880
9.1809E-03	1	11,1680

*Cell constant for sensor 1 (marked red)*



## References

- [1] RIEGER, Philip H. Electrochemistry. Englewood Cliffs, New Jersey 07632: Prentice-Hall, 1987. ISBN 0-13-248907-4.
- [2] Bard, Allen J.. "Electrochemical methods : fundamentals and applications / Allen J. Bard, Larry R. Faulkner." (1980).



Our visitor



Head Office – Strážek

### Contact us for more information

E-mail: [info@bvt.cz](mailto:info@bvt.cz)

Mob.: +420 775 605 538

Visit us on the webpage [www.bvt.cz](http://www.bvt.cz)

**BVT Technologies, a.s.**

**Company head office:**

Strážek 206, 592 53 Strážek, Czech Republic



[www.bvt.cz](http://www.bvt.cz)

[info@bvt.cz](mailto:info@bvt.cz)