

**SIMPLE CONNECTOR FOR ELECTROCHEMICAL SENSORS
 TYPE AC1, AC4, AC11, AC13, AP1, CC1, CC2, TS1**

Type: KA1.S.*

Description

The connector enables the use of the biosensor based on the substrates AC1, AC4, AC11, AC13, CC1, CC2 and TS1 in classical electrochemical arrangement with stirred vessel.



Physical Parameters

Dimensions:
 Weight: 28 gms
 Length: 39 mm
 Width: 9 mm
 Cable length 1.5 m

Types of termination

Model	Cable	Termination	Evaluating Units
KA1.S.S	Three core Shielded cable	Single conductors	any device
KA1.S.1	Three core Shielded cable	Banana plugs 4mm	any device
KA1.S.2	Three core Shielded cable	7 poles BVT connector	any device
KA1.S.3 TRIAD	Three core Shielded cable	Triad01 Palminstruments	Palmsens, EmStat2
KA1.S.3 LEMO4	Three core Shielded cable	LEMO 4 pins Palminstruments	Palmsens3, EIS, EmStat3, 3+
KA1.S.3 LEMO5	Three core Shielded cable	LEMO 5 pins Palminstruments	Palmsens4
KA1.S.4	Three core Shielded cable	BNC connectors	eDAQ
KA1.S.5	Three core Shielded cable	Banana plugs 2mm	any device

Ordering Information


- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA1.S connector is 4 weeks from receipt of order
- Delivery time for non-standard KA1.S connector depends on final technical specification of order

Example of Order


- 5 pieces - KA1.S.3 LEMO4

Internal Wiring

KA1.S.S



KA1.S.1



KA1.S.2




KA1.S.3 TRIAD, KA1.S.3 LEMO4
 KA1.S.3 LEMO5


TRIAD



LEMO



KA1.S.4



KA1.S.5

