


## SIMPLE CONNECTOR FOR ELECTROCHEMICAL SENSORS TYPE AC6

Type: KA6.S.\*

### Description

The connector enables the use of the biosensor based on the substrates AC6, in classical electrochemical arrangement with stirred vessel.



### Physical Parameters

#### *Dimensions:*

Weight: 35 gms  
 Length: 6 mm  
 Width: 8 mm  
 Cable length 1.5 m

### Types of termination

Model	Cable	Termination	Evaluating Units
KA6.S.S	Four core Shielded cable	Single conductors	any device
KA6.S.1	Four core Shielded cable	Banana plugs 4mm	any device
KA6.S.4	Four core Shielded cable	BNC connectors	any device
KA6.S.5	Four core Shielded cable	Banana plugs 2mm	any device

### Ordering Information


- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA6.S connector is 4 weeks from receipt of order
- Delivery time for non-standard KA6.S connector depends on final technical specification of order

### Example of Order


- 5 pieces - KA6.S.S

**Internal Wiring**


KA6.S.S



KA6.S.1



KA6.S.4



KA6.S.5

