

BVT pmw regulator for ST1

Type: BVT regulator

Description

- Universal BVT regulator for stirrer ST1.
- The input voltage 5 V
- USB port enables the use as USB device (PC, Power bank, Mobile charger.
- The output voltage 0-5V 2mm banana

Physical Parameters

Dimensions:

Weight: 55 gms
Length: 85 mm
Width: 59 mm
Height: 53 mm
USB cable length: 1m

Delivery ST1 contains

- BVT pmw regulator
- USB cable



Ordering information

- The order is specified by the whole product code
- Minimum order quantity - 1 BVT pmw regulator
- Delivery time for standard BVT pmw regulator is 4 weeks from receipt of order
- Delivery time for non-standard BVT pmw regulator depends on final technical specification of order
- Example of Order
 - 2 pieces BVT pmw regulator