

SIMPLE CONNECTOR FOR ELECTROCHEMICAL SENSORS TYPE AC8

Type: KA8.S.*


Description

The connector enables the use of the biosensor based on the substrates AC8 in classical electrochemical arrangement.

Physical Parameters

Dimensions:

Weight: 35 gms
 Length: 6 mm
 Width: 8 mm
 Cable length 1 m



Types of termination

Model	Cable	Termination	Evaluating Units
KA8.S.S	Four core Shielded cable	Single conductors	any device
KA8.S.1	Four core Shielded cable	Banana plugs 4 mm	any device
KA8.S.4	Four core Shielded cable	BNC connectors	any device
KA8.S.5	Four core Shielded cable	Banana plugs 2 mm	any device

Ordering Information


- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA8.S connector is 4 weeks from receipt of order
- Delivery time for non-standard KA8.S connector depends on final technical specification of order

Example of Order


- 5 pieces - KA8.S.5

Internal Wiring

KA8.S.S



KA8.S.1



KA8.S.4



KA8.S.5

