

CABLES FOR ELECTRODES

Type: CABLE.*.*


Description

Cables for BVT makes


Physical Parameters

Length: 100 cm

| Model | Cable | End1 | End2 | Evaluating Units |
|---------------|----------------------------|---------------------|----------------------|-----------------------------|
| CABLE.S.Red | One core Shielded cable | 2 mm banana plug | Single conductors | Conic electrode WCEc |
| CABLE.S.Black | One core Shielded cable | 2 mm banana plug | Single conductors | Conic electrode ACEc, |
| CABLE.S.Blue | One core Shielded cable | 2 mm banana plug | Single conductors | Conic electrodes RCEc |



| Model | Cable | End1 | End2 | Evaluating Units |
|---------------|-------------------------|------------------|-----------------|-----------------------|
| CABLE.1.Red | One core Shielded cable | 2 mm banana plug | 4mm banana plug | Conic electrode WCEc |
| CABLE.1.Black | One core Shielded cable | 2 mm banana plug | 4mm banana plug | Conic electrode ACEc, |
| CABLE.1.Blue | One core Shielded cable | 2 mm banana plug | 4mm banana plug | Conic electrodes RCEc |



| Model | Cable | End1 | End2 | Evaluating Units |
|--------------|---------------------------|-------------------|----------------------|---------------------|
| CABLE.2.BVT7 | Three core Shielded cable | 2 mm banana plugs | 7 poles BVT Conector | EmStat 2, PalmSens2 |



| Model | Cable | End1 | End2 | Evaluating Units |
|---------------|---------------------------|-------------------|-----------------|--|
| CABLE.3.TRIAD | Three core Shielded cable | 2 mm banana plugs | Triad conector | EmStat 2, PalmSens2 |
| CABLE.3.LEMO4 | Three core Shielded cable | 2 mm banana plugs | Lemo 4 conector | EmStat3, EmStat3 Blue, MultiEmstat3, PalmSens3 |
| CABLE.3.LEMO5 | Five core Shielded cable | 2 mm banana plugs | Lemo 5 conector | PalmSens4 |



| Model | Cable | End1 | End2 | Evaluating Units |
|----------------------|----------------------------|---------------------|------|-----------------------------|
| CABLE.4.Red | One core Shielded cable | 2 mm banana plug | BNC | Conic electrode WCEc |
| CABLE.5.Black | One core Shielded cable | 2 mm banana plug | BNC | Conic electrode ACEc, |
| CABLE.5.Blue | One core Shielded cable | 2 mm banana plug | BNC | Conic electrodes RCEc |



| Model | Cable | End1 | End2 | Evaluating Units |
|----------------------|----------------------------|---------------------|--------------------|-----------------------------|
| CABLE.5.Red | One core Shielded cable | 2 mm banana plug | 2mm banana plug | Conic electrode WCEc |
| CABLE.5.Black | One core Shielded cable | 2 mm banana plug | 2mm banana plug | Conic electrode ACEc, |
| CABLE.5.Blue | One core Shielded cable | 2 mm banana plug | 2mm banana plug | Conic electrodes RCEc |

