

## Mini RDE tips

---

Type: Mini RDE tips





### Description

BVT provides a wide range of RDE tips, composed from different materials from gold, platinum and copper to glassy carbon and working electrode diameters from 2 to 5 mm (as listed below).

The electrode material is inserted into a body of the Rotating disc electrode (RDE) to ensure excellent chemical compatibility and mechanical integrity.




The electrode tip materials for Mini RDE are defined by 1RDE.W1 2 mm, 2RDE.W\* 2.5 mm, 3RDE.W4 3 mm, 3RDE.W4 5 mm

<b>1RDE.W1 2 mm</b>		<b>2RDE.W* 2.5mm</b>	
<b>1 = Wire in glass</b>		<b>2 = Screen printed electrode</b>	
RDE = Rotating disc electrode		RDE = Rotating disc electrode	
<b>W = Work electrode</b>		<b>W = Work electrode</b>	
<b>1 - Platinum</b>		<b>1 - Platinum</b>	
Diameter of working electrode: 2 mm		<b>2 - Gold</b>	
		<b>3 - Copper</b>	
		<b>4 - Glassy carbon</b>	
		Diameter of working electrode: 2.5 mm	
<b>3RDE.W4 3 mm</b>		<b>3RDE.W4 5 mm</b>	
<b>3 = Tip in a peek casing</b>		<b>3 = Tip in a peek casing</b>	
RDE = Rotating disc electrode		RDE = Rotating disc electrode	
<b>W = Work electrode</b>		<b>W = Work electrode</b>	
<b>4 - Glassy carbon</b>		<b>4 - Glassy carbon</b>	
Diameter of working electrode: 3 mm		Diameter of working electrode: 5 mm	

### Device Advantages

- Easy replacement of electrode tips for RDE using the take out tool device
- The PEEK material of the Mini RDE tips ensures high chemical resistance



### Additional accessory equipment

Mini RDE	Mini rotating disc electrode
Mini RDE Kit	Mini rotating disc electrode Kit
TC4, TC5, TC6, TC7	Electrochemical glass cells
WCEc	External Working Classic Electrode Conic (+cables)
ACEc	External Auxiliary Classic Electrode Conic (+cables)
PalmSens4	PalmSens - potentiostat
AC1, AC2, AC4	Electrochemical sensors
KA1.C.S	Conic connector for electrochemical sensors
STK.S-C	The customer starter kit - set of twenty different electrochemical sensors to find the best one that suits your application

### Ordering/Price Information

- Delivery time is 8 weeks from receipt of order
- Example of order: 2 pieces 3RDE.W4 3 mm, 3 pieces 2RDE.W4 2.5mm