

SIMPLE CONNECTOR FOR ELECTROCHEMICAL SENSORS TYPE AC6

Type: KA6.S.*


Description

The connector enables the use of the electrochemical sensors and biosensors based on the substrates AC6, in classical electrochemical arrangement with stirred vessel.

Physical Parameters

Dimensions:

Weight: 29 gms
 Length: 30 mm
 Width: 14 mm
 Cable length 1 m



Types of termination

Model	Cable	Termination	Evaluating Units
KA6.S.S	Four core Shielded cable	Single conductors	any device
KA6.S.1	Four core Shielded cable	Banana plugs 4 mm	any device
KA6.S.4	Four core Shielded cable	BNC connectors	any device
KA6.S.5	Four core Shielded cable	Banana plugs 2 mm	any device

Ordering Information


- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA6.S connector is 4 weeks from receipt of order
- Delivery time for non-standard KA6.S connector depends on final technical specification of order

Example of Order


- 5 pieces - KA6.S.S

Internal Wiring


KA6.S.S




KA6.S.1



KA6.S.4



KA6.S.5



R - Right - AUXILIARY 1
 L - Left - AUXILIARY 2