

AUXILIARY CLASSIC ELECTRODE CONIC WITH PLATINUM PLATE

Type: ACEc.Pt plate

Description

ACEc is a glass tube with **platinum plate of guaranteed purity 99.99 %**. Standard connection is with a 2 mm banana plug. Auxiliary classic electrode with platinum plate (ACEc.Pt plate) is designed for electrochemical measurements.

The size of the platinum sheet can be adjusted at the customer's request.

Physical Parameters

*Dimensions of Pt plate: 15 x 5 mm Pt
(Pt thickness approx. 0.21 mm)*

Length: 61 mm
Diameter of glass tube: 6 mm
Weight: 19.5 gms

Electrode material

ACEc.Pt plate (Auxiliary Classical Electrode Conic with platinum plate) - Pt

Ordering Information

- The order is specified by whole product code
- Minimum order quantity - 1 electrode
- Delivery time for standard ACEc electrode is 4 weeks from receipt of order
- Delivery time for non-standard ACEc electrode depends on final technical specification of order

Example of Order: - 10 pieces - ACEc.Pt plate

Connecting cable

Standard connection is with a 2 mm banana plug.

For choosing another type of termination (connector), look at the cables catalog and choose the one you need.

