

DESICCATOR

Type: Desiccator

Description

An airtight container which maintains an atmosphere of low humidity through the use of a suitable drying agent which occupies the bottom part of the desiccator. It is for the storage of dry objects that must not be exposed to the moisture normally present in the atmosphere.

Physical parameters

Inner Dimensions:

Weight: 9000 gms
Height: 420 mm
Width: 302 mm
Length: 324 mm

Three shelves with 85 mm space.



Ordering information

- Minimum order quantity - 1 desiccator
- Delivery time for standard desiccator is 4 weeks from receipt of order
- Delivery time for a non-standard desiccator depends on the final technical specification

Example of Order

- 2 pieces - Desiccator

[Datasheet: Desiccator](#)